

Intermountain Visage Radiology Viewer Client Installer & First Use Guide - Windows

Note: Viewing radiology images within the Visage Radiology Client requires specific access to be present on your Intermountain account. If you had previous access with the IMPAX/ResMD clients, your account has been provisioned for Visage. If you cannot login, please initiate a support ticket using the Affiliated Support Email hco@intermountain.co.ihc.com Accessing images via the Visage Client requires a Secure Connection e.g. Hardware, Firewall or Software VPN (SecureAccess).


1. In addition to a Secure Connection, 2-Factor Authentication is required for Visage Viewer login. If you don't have it on your account, go to the Affiliated Support site ihhelp.info -> 2-Factor Auth / Portal / Email Resources and follow the Ping ID 2-Factor Auth Setup Guide.

2. At the same Affiliated Support page ihhelp.info under Windows User Resources go to Option 5: Visage Radiology Viewer and click, download and install the Visage Client Installer, as circled to the right.


Option 5

- Visage Radiology Viewer Client Software
- Visage Client Installer (REMOVES the IMPAX Client)
- Visage Installation and First-Use Configuration Guide
- Visage Saving Patient Lists Guide
- Visage Imaging YouTube Training Channel


3. Assuring that your Secure Connection (Hardware/Firewall/SecureAccess) is in place/connected, double-click the IHC Visage Radiology desktop shortcut, which opens a browser to login/authenticate. On the first launch, a browser window like the one below will appear. Check the Always allow box as circled and then click Open Visage 7.



4. Enter your Intermountain userid (username) and, if this is your personal device, check the Remember my username box then click Next. If this is your first time logging in on the device, enter your password then click Sign In. The Ping ID authentication window will appear where you will click Sign On again. Once you authenticate on your smartphone, it will log into Visage.



5. To avoid this message every time you open Visage, you will want to check Private networks then Allow Access, as shown below:



6. When you see this screen, you have successfully launched the Visage Viewer:

The screenshot displays the Visage Viewer application interface. At the top, there is a menu bar with 'File', 'View', and 'Help' options, and a user profile indicator 'userid@imp01'. Below the menu bar, there are icons for 'View', 'Add to View', and 'Retrieve'. A '+ New Preset' button is visible on the left side.

The main interface is divided into several sections:

- Query Section:** Contains a search bar and a 'Query' button.
- Patient Information Section:** Includes fields for 'Name' (with a sub-field for 'Surname, Firstname'), 'Date of Birth' (format MM/DD/YYYY), and 'ID'.
- Filtering Section:** Includes a 'Type' dropdown set to 'Study Date', 'From' and 'To' date pickers (showing '2021-Oct-21 00:00' and '2021-Nov-05 23:59' respectively), and buttons for 'Today', 'Yesterday', '1 Week', and '2 Weeks'.
- Modality Section:** A checked 'Modality' section with radio buttons for 'CT', 'MR', 'PT', 'US', 'CR', 'DX', 'MG', and 'XA', and an 'Others' text input field.
- Fields Section:** Includes dropdown menus for 'Accession Number', 'Institution Name', and 'Study Labels', and a 'Study States...' button.
- Table Section:** A table titled 'My Recent Sessions' with columns for 'Patient Name', 'Session Name', 'Creator', 'Creation Time', and 'Studies'. The table is currently empty.
- Buttons:** A 'Refresh' button is located below the table, and a 'Show My Recent Sessions' button is in the bottom right corner.